

schülke +

mikrozid®


Surface disinfection against COVID-19,
when compromise is not an option

Effective against
coronaviruses
including
SARS-CoV-2



we protect lives
worldwide

SARS-CoV-2 is the virus that causes COVID-19



Contaminated hands

*SARS-CoV-2 can remain
on human skin for
around 9 hours*

[Hirose, 2020]



Transmission of SARS-CoV-2

*Spreads between people, when an infected person coughs,
sneezes or talks [WHO, 2020].*



Contaminated surfaces

SARS-CoV-2 found to be stable on:


Plastic and stainless steel = up to 72 hours

Cardboard = up to 24 hours

Copper = up to 4 hours

[van Doremalen, 2020]

*Human coronaviruses can remain infectious
on inanimate surfaces at room temperature
for up to 9 days [Kampf 2020]*



'Respiratory secretions or droplets expelled by infected individuals can contaminate surfaces and objects, creating fomites (contaminated surfaces). Viable SARS-CoV-2 can be found on those surfaces for periods ranging from hours to daysin particular at high concentration in health care facilities where COVID-19 patients were being treated.'

[WHO, 2020]

Cleaning without compromise

Public Health England recommends:

'The frequency of cleaning across all risk pathways should be increased during the pandemic to at least twice daily. Frequently touched sites/points, in particular waiting rooms/dental chair, should be cleaned between patients. A....disinfectant effective against viruses, bacteria and fungi to EN standard 14476 for viricidal activity should be used.' [PHE, 2021]



References

- Hirose R et al, Survival of SARS-CoV-2 and influenza virus on the human skin. Importance of hand hygiene in COVID-19, Clinical Infectious Diseases, 2020.
- Public Health England, COVID-19: infection prevention and control dental appendix, updated 21st January 2021
- van Doremalen N, Bushmaker T, Morris DH, et al. Aerosol and surface stability of SARS-CoV-2 as compared with SARS-CoV-1. *N Engl J Med.* 2020 Apr 16;382(16):1564-7
- Kampf G, Todt D, Pfaender S, Steinmann E. Persistence of coronaviruses on inanimate surfaces and their inactivation with biocidal agents. *Journal of Hospital Infection.* VOL 104, ISSUE 3, P246-251, 2020
- World Health Organization. Transmission of SARS-CoV-2: implications for infection prevention precautions – scientific brief. 2020 [internet publication] <https://www.who.int/news-room/commentaries/detail/transmission-of-sars-cov-2-implications-for-infection-prevention-precautions>

Selecting the ideal disinfectant [Rutala, 2014]



References

- Rutala WA, Weber DJ. Selection of the ideal disinfectant. *Infect Control Hosp Epidemiol.* 2014 Jul;35(7):855-65. doi: 10.1086/676877. Epub 2014 May 29. PMID: 24915214.

mikrozid®



limited	
VIR	30s
BAC	60s
YEAST	60s
TB	60s
Alcohol	

mikrozid® alcohol free



limited	
VIR	60s
BAC	60s
YEAST	60s
TB	-
Alcohol-free	

Material Compatibility

Rubber
Plastics
Metals

Material Compatibility


Rubber
Plastics
Metals
Leather
Plexi-glass

mikrozid® universal



limited	
VIR	15s
BAC	60s
YEAST	60s
TB	5 mins

mikrozid® products
are virucidal against
enveloped viruses
including coronaviruses
and SARS-CoV-2



Material Compatibility	
	Rubber Plastics Metals Leather Plexi-glass

Product Features		Ready To Use Wipes	mikroziid® universal	mikroziid® alcohol free
Composition	Active Ingredients			
		100g agent solution contains: - 25 g Ethanol (94 %) - 35 g Propan-1-ol	100g agent solution contains: - 12.6 g Ethanol (94 %) - 17.4 g Propan-2-ol - < 5% anionic surfactants	100g agent solution contains: quaternary ammonium compounds: - 0.26 g Alkyl(C12-C16) dimethylbenzyl ammonium chloride - 0.26 g Didecyldimethylammonium chloride - 0.26 g Alkyl(C12-C14) ethylbenzyl ammonium chloride
Contains alcohol?		✓	✓	✗
Metals (e.g. Stainless steel, copper, brass)		✓	✓	✓
Rubbers (e.g. Natural rubber, silicone)*		✓	✓	✓
Material Compatibility	Plastics (e.g. PE, PET, PS, PC)*	✓	✓	✓
Plexi-glass (e.g Acrylic)		✗	✓	✓
Leather		✗	✓	✓
Bactericidal		1 min	1 min	Clean condition: 1 min Dirty condition: 2 mins
Yeasticidal (e.g <i>Candida albicans</i> , <i>Candida auris</i>)		1 min	1 min	Clean condition: 1 min Dirty condition: 2 mins
Fungicidal (e.g. <i>Aspergillus</i> species)		clean: 2 mins	-	-
Mycobactericidal (e.g. <i>Mycobacterium terrae</i>)		1 min	5 mins	-
Antimicrobial Spectrum & ContactTime*	Enveloped virus (e.g. HIV, HBV, HCV, Coronaviruses)	30 sec	15 sec	1 min
	Non-enveloped virus (Rotavirus)	n/a	15 sec	1 min
	Non-enveloped virus (Adenovirus)	30 sec	-	-
	Non-enveloped virus (Norovirus)	1 min	clean: 30 sec; dirty: 1 min	15 mins
Surface Coverage		0.5m ² (Jumbo Tub)	1m ² (Premium Pack)	0.5m ² (Jumbo Tub)

*For more information, please contact the relevant salesperson or email to mail.uk@mail.schuelke.com

Legend: ✓ - Yes, ✗ - No

How to use

Jumbo Tub



Wash hands thoroughly.



Put on gloves.



Invert the tub 3 times.



Open lid and pull out the wipes.




Wipe the surface with an "S" pattern.



Ensure complete wetting of the surface throughout the exposure time.



Dispose the used wipes.



Make sure that the lid is securely closed.

Liquid



Wash hands thoroughly.




Put on gloves.



Spray liquid into wipe



Wipe the surface with an "S" pattern.



Alternatively spray liquid directly on to uneven surfaces




Ensure complete wetting of the surface throughout the exposure time.



Dispose the used wipes.

Premium Pack



Wash hands thoroughly.



Put on gloves.



Open pack and pull out the wipes.




Wipe the surface with an "S" pattern.




Ensure complete wetting of the surface throughout the exposure time.



Dispose the used wipes.



Make sure that the lid is securely closed.



Make sure to store pack on its side

Product Summary

The plus of schülke – mikrozid® range ready-to-use wipes

- tested to EN 14476 standards = efficacy against enveloped viruses including coronaviruses.
 - high quality non-woven material.
 - certified to EN 16615 the highest level of testing for antimicrobial wipes.
- → →
- effective in high-risk and high-touch areas for COVID-19.
 - tear-resistant and holds optimum disinfectant for effective cleaning.
 - the test examines the efficacy of the wipe as a whole i.e. the wipe plus the disinfectant component.



mikrozid®

limited
VIR 30s

Unique alcohol formulation for virucidal protection in the cleaning and disinfection of hard surfaces

BAC 60s

- blend of two alcohols for a very high performance against viruses within one minute contact time
- cleaning and disinfection, two complete actions in one product
- added surfactants for boosted cleaning performance
- effective against bacteria (inc. TB), yeast and viruses (enveloped viruses inc. Coronaviruses, HBV, HIV, HCV)
- quick drying and residue free



Pack size	Article n°
200 jumbo wipes tub	109 214
200 jumbo wipes refill	7000 1450
150 wipes tub	109 216
150 wipes refill	70002065
1 L bottle	109210
5 L canister	109211




mikrozid® alcohol free

limited
VIR 60s

Alcohol free, dual action cleaner and disinfectant for alcohol sensitive surfaces

BAC 60s

- ideal for all surfaces including including alcohol sensitive materials
- added surfactants for boosted cleaning performance
- cleaning and disinfection, two complete actions with one product
- effective against bacteria, yeast and viruses (enveloped viruses inc. Coronaviruses, HBV, HIV, HCV)



Pack size	Article n°
200 jumbo wipes tub	165 723
200 jumbo wipes refill	7000 1449
1 L bottle	165 720
5 L canister	165 721



mikrozid® universal

limited
VIR 15s


Uniquely versatile, low-alcohol formulation for the rapid cleaning and disinfection of surfaces & non-invasive medical devices

BAC 60s

- unique low alcohol formulation - blend of two alcohols
- dual action - cleaning and disinfection with one product
- added surfactants for outstanding cleaning
- compatible with a wide range of materials including: leatherette, tablets and displays
- effective against bacteria (inc. TB), yeast and viruses (enveloped viruses inc. Coronaviruses, HBV, HIV, HCV)
- dermatologically tested



Pack size	Article n°
100 wipes	7000 0038
1 L liquid	7000 0998



schülke worldwide:

Australia

Schülke Australia Pty Ltd
Macquarie Park NSW 2113
Phone +61 2 8875 9300
Fax +61 2 8875 9301

Austria

Schülke & Mayr Ges.m.b.H.
1070 Vienna
Phone +43 1 523 25 01 0
Fax +43 1 523 25 01 60

Brazil

Vic Pharma Indústria e
Comércio Ltda.
Taquaratinga/SP – CEP
15900-000
Phone +16 3253 8100
Fax +16 3253 8101

Czech Republic

Schülke CZ, s.r.o.
73581 Bohumín
Phone +420 558 320 260
Fax +420 558 320 261

France

Schülke France SARL
92250 La Garenne-
Colombes
Phone +33 1 42 91 42 42
Fax +33 1 42 91 42 88

India

Schülke India Pvt. Ltd.
New Delhi 110044
Phone +91 11 30796000
Fax +91 11 42595051

Japan

Schülke Japan LLC
Shinagawa East One Tower
2-16-1 Konan, Minato-ku,
Tokyo 108-0075
Phone +81 3 6894 3335
Fax +81 3 3453 0207

Malaysia

Schülke & Mayr (Asia) Sdn Bhd.
47301 Petaling Jaya, Selangor
Phone +60 3 78 85 80 20
Fax +60 3 78 85 80 21

Netherlands

Schülke & Mayr Benelux B.V.
2031 CC Haarlem
Phone +31 23 535 26 34
Fax +31 23 536 79 70

New Zealand

Schülke New Zealand Limited
Auckland 1010
Phone +61 2 8875 9300
Fax +61 2 8875 9301

Poland

Schülke Polska Sp. z o.o.
02-305 Warszawa
Phone +48 22 11 60 700
Fax +48 22 11 60 701

Singapore

Schülke & Mayr (Asia) Pte. Ltd.
Singapore 159410
Phone +65 62 57 23 88
Fax +65 62 57 93 88

Slovakia

Schülke SK, s.r.o.
97101 Prievidza
Phone +421 46 549 45 87
Fax +420 558 320 261

Switzerland

Schülke & Mayr AG
8003 Zurich
Phone +41 44 466 55 44
Fax +41 44 466 55 33

United Kingdom

Schülke & Mayr UK Ltd.
Sheffield S9 1AT
Phone +44 114 254 35 00
Fax +44 114 254 35 01

... plus our international distributors

schülke -

Schülke & Mayr UK Ltd.

Sheffield S9 1AT | United Kingdom
Phone +44 114 254 35 00
Fax +44 114 254 35 01
www.schuelke.com

 www.youtube.com/schuelkeChannel

 www.facebook.com/myschuelke

 twitter.com/schuelkeUK